PMMA INTRAOCULAR LENS

DEVICE DESCRIPTION

Intraocular lenses are optical implants for the replacement of the crystalline lens for the visual correction of aphakia. All lenses are lathe cut and tumble polished from the finest polymethylmethacrylate (PMMA) with ultraviolet (UV) absorber.

The material has a refractive index of 1.49, a specific gravity of 1.19 and a transmission of visible light of > 91 %. The cut- off wavelength of 10 % transmission for UV light is 372 nm for a +4 diopter lens and 388 nm for a +34 diopter lens.

All lenses are supplied sterile and the inner package is terminally sterilized and should be opened under sterile conditions.

MODE OF ACTION / INDICATIONS

Posterior IOLs are a refracting medium to replace the natural crystalline lens in correcting aphakia where the cataractous lens has been removed by extracapsular extraction methods.

CAUTION

Patients with the following conditions may not be suitable candidates for IOL due to interference with or exacerbation of pre-existing conditions:

- Congenital Bilateral Cataracts
- Recurrent Anterior / Posterior Inflammation
- Microphthalmos
- Chronic Angle-Closure Glaucoma
- Persistent/Significant bleeding or vitreous loss
- Monocular vision
- Glaucoma (Non-controllable)
- Corneal Endothelial Dystrophy
- Proliferative diabetic Retinopathy
- Anyone under 18 years of age

WARNINGS

All surgical procedures carry certain risks. IOL implantations carry the following potential risks and complications. Lens discoloration, non pigment precipitates, corneal endothelial damage, high IOP, infection, retinal detachment vitritis, cystoid macular edema, corneal edema, papillary membranes, flat anterior chamber, iris prolapsed, phpopoyon, and transient glaucoma.

- The safety of IOL implantation has not been substantiated in patients with any pre-existing ocular medical conditions.
- Long term effects of IOL implants have not been determined. Use caution with all IOLs.
- Some patients with macular degeneration may NOT have improved visual outcomes.
- Patients with pre-existing conditions may not get as good of results as those without.
- Patients should always be carefully monitored postoperatively, especially those with any post-operative complications.
- Posterior IOLs should not usually be implanted in the Anterior chamber of the patient.

PRECAUTIONS

- Do NOT resterilize IOLs
- Do NOT store in direct sunlight
- Use only BSS or equivalent

When reused, there is an increased RISK for material damage and increased risk of infection.

For Single Use Only!

REPORTING

All adverse reactions that may be lens related should be reported to the manufacturer in order to document long term effects of IOLs. If applicable contact at:

i-Medical® Ophthalmic International Heidelberg GmbH

Markircher Straße 7 · D-68229 Mannheim · Germany
Phone: +49 (0) 621-4844 90-0 · Fax: +49 (0) 621-4844 90-20
www.imedical.de

LENS POWER CALCULATION

The physician should determine pre-operatively the power of the IOL to be implanted. IOL calculation methods are in the references at the end of this insert.

DIRECTIONS

Remove the lens tray from the pouch by peeling open the seal in a sterile environment. Open the case by either turning the cape counter-clockwise or sliding back the cover, to present the lens for forceps extraction.

NOTE: Lenses are plastic and may carry a static charge the will attract dust and particulates. Remove such particles by rinsing in BSS equivalent prior to implantation.

Remove the lens by grasping the optic on the opposite sides, taking care not to scratch the optical surface. Locking forceps or needle holders should never be used.

PMMA IOL haptics are designed to compress with the plane of the lens ONLY. Do not attempt to expand the haptics or flex them out of the plane of the optic and do not torque or twist the lens.

Prior to implantation examine the lens.

IMPLANT CARD

A patient implant identification card is provided with each IOL. Patients should be instructed to keep the card as a permanent record of the implant and to give it to their eye care practitioner.

EXPIRATION DATE

The expiration date on the package is the sterility expiration date.

HOW SUPPLIED

i-Medical® lenses are supplied sterile in a dry, heat/sealed package that is terminally sterilized. Open in sterile environment only.

BIBLIOGRAPHY

The following list of articles, books and publication is only a small fraction of material available on IOLs and implantation. Physicians should carefully study, complete courses and feel comfortable with IOLs prior to their use with patients.

GENERAL

Apple, D.J.; Evolution of Intraocular Lenses, Ophthalmology, 1985, 1-23

Holiday, J.T.; The Ideal Intraocular lens, CLAO Journal, 1983, 15-19.

Jaffee, N.S.; The way things were and are, changing indications for IOL Implantation; Ophthalmology, 1983, 381.

Shepard, D.D.; The IOL lens Manual, Shepard Inc., 1982.

Mc Donell, P.J.; Green, W.R.; Pathology of IOLs in 33 Eyes Examined Postmortum, Ophthalmology, 1983, 386.

Ridley, Harold; The Origin and Objectives of IOL lenticular

Implants; American Academy of Ophthalmology and Otolaryngology Transactions; 1976,65.

POSTERIOR CHAMBER IOLS

Gills, J.P.; My method of Extracapsular Cataract Extraction with Implantation of a Posterior Chamber IOL, Ophthalmic Surgery, 1985, 386-392.

Harris, W.W.; Kogan, I.; Retrospective Study of Retina Complications with Posterior Chamber IOLs, Contact/IOL Medical Journal, 1981, 341-344.

Jaffe, N.S.; Current Concepts in Posterior Chamber lens Technology, American Implant Society Journal, 1985, 456.

Pearce, J.L.; Sixteen months experience with 140 posterior chamber IOLs, British Journal of Ophthalmology, 1977, 310.

Binkhorst, C.D.; Manshot, W.A.; The Mechanism of Capsular Fixation of IOLs in the Human Eye, American Implant Society Journal, 1976, 110.

Alpar, J.J. Guidelines for Indication of IOL Implantation and the Selection of the Proper Lens, Ophthalmic Surgery, 1995.

POWER CALCULATIONS

Binkhorst, R.D., IOL Power Calculation Manual, New York, Richard D. Binkhorst, 1978.


Donzis, P.B., An IOL lens Formula for Short, Normal and Long Eyes, CLAO Journal, 1985, 95.

Gills, J.P.; Minimizing Postoperative Refractive Error, Contact Lens, 1980, 56.

Holladay, J.T.; Improving the Predictability of IOL Power Calculations, Arch Ophthalmology, 1986, vol. 104, pg. 538.

Retzlaff, J.A.; Manual of Implant Power Calculation, SRK Formula, 1981.

Richards, S.C.; Clinical Evaluation of Six IOL Calculation Formulas, 1985.

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